


... INSTRUMENTATION



INSTRUMENTS DE CHIRURGIE | INSTRUMENTS ROTATIFS CHIRURGICAUX



FRAISES CHIRURGICALES TUNGSTÈNE

■ ■ ■ H141



T00517	ISO 500 104 001 003 010 H141 104 010	010	-	-	-	-	-	-	-	-	
T00518	ISO 500 104 001 003 014 H141 104 014	-	014	-	-	-	-	-	-	-	
T00519	ISO 500 104 001 003 018 H141 104 018	-	-	018	-	-	-	-	-	-	
T00520	ISO 500 104 001 003 023 H141 104 023	-	-	-	023	-	-	-	-	-	
T00521	ISO 500 104 001 003 027 H141 104 027	-	-	-	-	027	-	-	-	-	
T00522	ISO 500 104 001 003 031 H141 104 031	-	-	-	-	-	031	-	-	-	
T00523	ISO 500 104 001 003 035 H141 104 035	-	-	-	-	-	-	035	-	-	
C03027	ISO 500 104 001 003 040 H141 104 040	-	-	-	-	-	-	-	040	-	
C03028	ISO 500 104 001 003 050 H141 104 050	-	-	-	-	-	-	-	-	050	
T00525	ISO 500 105 001 003 023 H141 105 023	-	-	-	023	-	-	-	-	-	
T00526	ISO 500 105 001 003 027 H141 105 027	-	-	-	-	027	-	-	-	-	
T00527	ISO 500 105 001 003 031 H141 105 031	-	-	-	-	-	031	-	-	-	
T00528	ISO 500 105 001 003 035 H141 105 035	-	-	-	-	-	-	035	-	-	
C03029	ISO 500 105 001 003 040 H141 105 040	-	-	-	-	-	-	-	040	-	
C03030	ISO 500 105 001 003 050 H141 105 050	-	-	-	-	-	-	-	-	050	

■ ■ ■ COFFRET CHIRURGIE



C11010 **Assortiment :**
H151 314 016 - H151XL 315 016 - H162A 314 016 - H253 316 006


FG 315-316-317 : < 160.000 min.⁻¹ CA 204 : 5.000 min.⁻¹ à 50.000 min.⁻¹ PM 104-105-106 : 5.000 min.⁻¹ à 50.000 min.⁻¹




■ ■ ■ 33IL (FRAISE IMPLANTAIRE)



T10399	FG 316	ISO 500 316 415 007 010 HM33IL 316 010	-	-	-	-	-	-	-	010
--------	--------	--	---	---	---	---	---	---	---	-----

■ ■ ■ H1S



T00501	ISO 500 204 001 003 010 H1S 204 010	010	-	-	-	-	-	-	-	-	
T00502	ISO 500 204 001 003 012 H1S 204 012	-	012	-	-	-	-	-	-	-	
T00503	ISO 500 204 001 003 014 H1S 204 014	-	-	014	-	-	-	-	-	-	
T00504	ISO 500 204 001 003 016 H1S 204 016	-	-	-	016	-	-	-	-	-	
T00505	ISO 500 204 001 003 018 H1S 204 018	-	-	-	-	018	-	-	-	-	
T00506	ISO 500 204 001 003 021 H1S 204 021	-	-	-	-	-	021	-	-	-	
T00507	ISO 500 204 001 003 023 H1S 204 023	-	-	-	-	-	-	023	-	-	
T00509	ISO 500 205 001 003 018 H1S 205 018	-	-	-	-	018	-	-	-	-	
T00510	ISO 500 205 001 003 023 H1S 205 023	-	-	-	-	-	-	023	-	-	
T00511	ISO 500 205 001 003 027 H1S 205 027	-	-	-	-	-	-	-	027	-	
T00513	ISO 500 206 001 003 018 H1S 206 018	-	-	-	-	018	-	-	-	-	
T00514	ISO 500 206 001 003 023 H1S 206 023	-	-	-	-	-	-	023	-	-	
T00515	ISO 500 206 001 003 027 H1S 206 027	-	-	-	-	-	-	-	027	-	

COMMANDES : N°AZUR 0 810 001 801 MAIL chirurgie@dexter.fr





FRAISES CHIRURGICALES TUNGSTÈNE

■ ■ ■ H151



T00176	FG 314	ISO 500 314 199 336 016 H151 314 016	-	-	-	-	016
T00177	FG 316	ISO 500 316 199 336 016 H151XL 316 016	-	-	-	-	016
T00534	CA 204	Séparation de racine H151XL 204 ...	-	-	-	-	016
T00545	PM 104	ISO 500 104 199 336 016 H151 104 016	-	-	-	-	016
C03023	PM 106	ISO 500 106 199 336 016 H151XL 106 016	-	-	-	-	016

■ ■ ■ H157



T00181	FG 314	H157 314 012	-	-	-	-	012
--------	--------	---------------------	---	---	---	---	-----

■ ■ ■ H162A



T00188	FG 315	ISO 500 315 408 298 016 H162A 315 016	-	-	-	-	016
T05451	CA 205	ISO 500 205 408 298 016 H162A 205 016	-	-	-	-	016
T05450	PM 104	ISO 500 104 408 298 016 H162A 104 016	-	-	-	-	016

■ ■ ■ H166A



C03035	CA 206	ISO 500 206 409 298 021 H166A 206 021	-	-	-	-	021
C03034	PM 104	ISO 500 104 409 298 021 H166A 104 021	-	-	-	-	021

⌚ FG 315-316-317 : < 160.000 min.⁻¹ CA 204 : 5.000 min.⁻¹ à 50.000 min.⁻¹ PM 104-105-106 : 5.000 min.⁻¹ à 50.000 min.⁻¹

■ ■ ■ H161 LINDEMANN



T00186	FG 315	ISO 500 315 408 295 016 H161 315 016	-	-	-	-	016
T00541	PM 104	ISO 500 104 408 295 016 H161 104 016	-	-	-	-	016

■ ■ ■ H162 LINDEMANN



T00187	FG 315	ISO 500 315 408 297 016 H162 315 016	-	-	-	-	016
T00529	CA 205	ISO 500 205 408 297 016 H162 205 016	-	-	-	-	016
T00531	CA 206	ISO 500 206 408 297 016 H162 206 016	-	-	-	-	016
T00542	PM 104	ISO 500 104 408 297 016 H162 104 016	-	-	-	-	016

■ ■ ■ H166 LINDEMANN



T00543	PM 104	ISO 500 104 409 297 02 H166 104 021	-	-	-	-	021
--------	--------	---	---	---	---	---	-----

■ ■ ■ H253



T00191	FG 314	ISO 500 314 415 296 006 H253 314 006	-	-	-	-	006
--------	--------	--	---	---	---	---	-----




... INSTRUMENTATION

INSTRUMENTS DE CHIRURGIE | INSTRUMENTS ROTATIFS CHIRURGICAUX

FRAISES CHIRURGICALES TUNGSTÈNE


■ ■ ■ H254



T00192	 FG 314	ISO 500 314 415 296 010 H254 314 010	-	-	-	-	-	010
--------	--	--	---	---	---	---	---	-----

■ ■ ■ H269




T00194	 PM 104	ISO 500 104 199 295 016 H269 104 016	-	-	-	-	-	016
--------	--	--	---	---	---	---	---	-----

FRAISES CHIRURGICALES INOX

■ ■ ■ 82



GN4291	 PM 106	ISO 330 106 155 172 050 82 106 050	-	-	-	-	-	050
--------	--	--	---	---	---	---	---	-----




■ ■ ■ 85



GN4321	 PM 106	ISO 330 106 260 172 050 85 106 050	-	-	-	-	-	050
--------	--	--	---	---	---	---	---	-----




■ ■ ■ LINDEMANN



T00540	ISO 330 104 406 297 014 RF163 104 014	014	-	-	-	-	 PM 104
T10540	ISO 330 104 406 297 018 RF164 104 018	-	018	-	-	-	-
T20540	ISO 330 104 406 297 023 RF165 104 023	-	-	023	-	-	-
T30540	ISO 330 104 406 297 021 RF166 104 021	-	-	-	021	-	-
T40540	ISO 330 104 406 297 023 RF167 104 023	-	-	-	-	023	-
T50540	ISO 330 105 406 297 023 RF168 105 023	 PM 105	-	-	-	-	023
T60540	ISO 330 106 406 297 023 RF169 106 023	 PM 106	-	-	-	-	023

■ ■ ■ DINGMANN



T00549		= 1,5	4453 106 020	020
T10549		= 2,0	4455 106 028	028
T20549		= 1,5	4457 106 028	028

