




FRAISES CHIRURGICALES DIAMANTÉES

TREPANS

Ø < 045 : < 4.000 min.⁻¹ Ø > 045 : < 2.000 min.⁻¹

242 DIAMANTÉE

◆ 524 - Grain Moyen



C70242	ISO 806 106 001 524 018 DIAM 242 106 018	018	-	-	-	-	-	-	PM 106
T30242	ISO 806 106 001 524 023 DIAM 242 106 023	-	023	-	-	-	-	-	-
T40242	ISO 806 106 001 524 027 DIAM 242 106 027	-	-	027	-	-	-	-	-
T00242	ISO 806 106 001 524 030 DIAM 242 106 030	-	-	-	030	-	-	-	-
T10242	ISO 806 106 001 524 035 DIAM 242 106 035	-	-	-	-	035	-	-	-
T20242	ISO 806 106 001 524 040 DIAM 242 106 040	-	-	-	-	-	040	-	-
T50242	ISO 806 106 001 524 045 DIAM 242 106 045	-	-	-	-	-	-	045	-
T60242	ISO 806 106 001 524 050 DIAM 242 106 050	-	-	-	-	-	-	-	050

Os compact : ⚙️ FG 314 : 20.000 min.⁻¹ à 120.000 min.⁻¹

Os spongieux : ⚙️ FG 314 : 10.000 min.⁻¹ à 30.000 min.⁻¹

G1



T00196	FG 314	ISO 806 314 411 534 G1 314 021	-	-	-	-	-	-	021
---------------	--------	--	---	---	---	---	---	---	-----


C249



⚙️ FG 314 : 300.000 min.⁻¹ à 500.000 min.⁻¹ SANS SPRAY

C00199	FG 314	C249 314 016	-	-	-	-	-	-	016
C10199		C249 314 016	-	-	-	-	-	-	016

229




C00569	ISO 330 205 486 001 229 205 020	⊥ 3,0 020	-	-	-	-	-	CA 205
C00570	ISO 330 205 486 001 229 205 035	-	⊥ 4,5 035	-	-	-	-	-
C10570	ISO 330 205 486 001 229 205 040	-	-	⊥ 5,0 040	-	-	-	-
C40570	ISO 330 205 486 001 229 205 055	-	-	-	⊥ 6,5 055	-	-	-
C50570	ISO 330 205 486 001 229 205 060	-	-	-	-	⊥ 7,0 060	-	-
C00560	ISO 330 205 486 001 229 205 070	-	-	-	-	-	⊥ 8,0 070	-

Longueur = 10 mm

⊥ Diamètre extérieur (mm)

229L



C30580	ISO 330 205 555 001 229L 205 050	⊥ 6,0 050	-	-	-	-	-	CA 205
C40580	ISO 330 205 555 001 229L 205 055	-	-	-	-	-	⊥ 6,5 055	-

Longueur = 14 mm

⊥ Diamètre extérieur (mm)

229XL



C30590	ISO 330 205 555 001 229XL 205 050	⊥ 6,0 050	-	-	-	-	-	CA 205
---------------	---	--------------	---	---	---	---	---	--------

Longueur = 18 mm

⊥ Diamètre extérieur (mm)

